



RapidScan ST5-W Lateral Flow Assay Reader

Validated for multiple food and environmental lateral flow tests

Food or feed manufacturers, grain trading companies and laboratories using lateral flow devices as a rapid test method require to obtain quick and reliable results for contamination testing.

RapidScan ST5-W Lateral Flow Assay Reader is a highly flexible, easy to use and accurate imaging platform used for lateral flow assay analysis.

The instrument was designed specifically for field and in-process testing applications requiring qualitative to quantitative test results.

RapidScan ST5-W is unique to the market with its flexibility, multiple target analysis and its **compatibility with Eurofins Technologies' allergen, GMO, glyphosate and algal toxins lateral flow tests**, providing an integrated solution to customers.

Key Benefits

- User friendly: no training needed, easy to use
- Semi-quantitative results
- Portable for field testing
- Flexible: multiple target analysis
- Data storage capacity of up to 5000 test results
- Report exporting capability

RapidScan ST5-W Lateral Flow Assay Reader



- **Ease-of-use:** the reader is easy to set up and operate, applying an intuitive user interface with color touchscreen display.
- **Quantitative results** are obtained in a few minutes. Data can be stored on the reader then exported through a USB port.
- **Applications:** can be used with a wide range of allergen, GMO and glyphosate strip tests.



Allergen testing with RapidScan ST5-W

RapidScan is unique to the market for multiple allergen testing, compatible with Eurofins Technologies' growing range of allergen lateral flow devices. These rapid tests, with **high sensitivity and specificity**, are provided with all the necessary material to perform 20 individual analysis per kit, and validated on **food matrix, swabbing and rinse water**.

Article No.	Product Name
HU0030110	SENSIStrip Almond
HU0030111	SENSIStrip Beta-Lactoglobulin
HU0030112	SENSIStrip Casein
HU0030113	SENSIStrip Cashew
HU0030114	SENSIStrip Coconut
HU0030115	SENSIStrip Crustaceans (Tropomyosin)
HU0030116	SENSIStrip Egg White
HU0030117	SENSIStrip Fish
HU0030118	SENSIStrip Gluten
HU0030119	SENSIStrip Hazelnut
HU0030120	SENSIStrip Lupin
HU0030121	SENSIStrip Lysozyme

Article No.	Product Name
HU0030122	SENSIStrip Macadamia nut
HU0030123	SENSIStrip Milk Protein
HU0030124	SENSIStrip Molluscs (Tropomyosin)
HU0030125	SENSIStrip Mustard
HU0030126	SENSIStrip Ovalbumin
HU0030127	SENSIStrip Brazil nut
HU0030128	SENSIStrip Peanut
HU0030129	SENSIStrip Pecan nut
HU0030130	SENSIStrip Pistachio
HU0030131	SENSIStrip Sesame
HU0030132	SENSIStrip Total Soy
HU0030133	SENSIStrip Walnut

Please note that some LFD kits will only be available later in 2022.





GMO and glyphosate testing with RapidScan ST5-W

GMO rapid testing is especially important for **grain traders** and elevators, as a primary screening tool of grain lots at reception points.

Food manufacturers can perform GMO screening of raw ingredients like corn, soybean and canola before it is processed into a final product. **Seed testing laboratories** providing testing services for bulk grain analysis can quickly and easily identify the amount of GM contamination.

Our most requested kits for GMO rapid testing are available on the **RapidScan ST5-W** reader:

Article No.	Product Name	Product Description
EAID 101	CP4 EPSPS (RR) Quantitative LFS Kit Soybean & Canola	RapidScan (0.1%) Bulk Grain Soybean & Canola, 100 strips Sensitivity: 0.1%
EAID 102	CP4 EPSPS (RR) Quantitative LFS Kit Corn	Bulk Grain Corn, 100 strips Sensitivity: 0.1%

Further to GMO testing, food companies and grain processors can run **pesticide (glyphosate)** strip tests on RapidScan, to prevent that contaminated raw materials are being used in finished food products.

Article No.	Product Name	Testing Range
500095	Glyphosate, Dipstick, 20 tests	0.5 - 20 ppb (water)



Algal toxins testing with RapidScan ST5-W

For **public water systems, lake managers, government agencies, regional testing programs, communities and individual researchers** requiring near real time answers on water safety to protect public health against the threat of algal toxins, our fresh and marine water toxin dipstick tests, applicable on RapidScan, provide the best solution.

Article No.	Product Name	Testing Range
520017 520046	Microcystins, Dipstick, Finished Drinking Water, 20 tests	0 - 5 ppb 0.3 - 5 ppb
520020	Microcystins, Dipstick, Source Drinking Water with QuikLyse® Feature, 20 tests	0 - 5 ppb
520022	Microcystins, Dipstick, Recreational Fresh Water with QuikLyse® Feature, (EPA ETV), 20	0 - 10 ppb
520030	Cylindrospermopsin, Dipstick, with QuikLyse® Feature, 20 tests	0 - 10 ppb
520043	Anatoxin-a (VFDF), Dipstick, 20 tests	0 - 2.5 ppb
520045	Saxitoxin (PSP), Dipstick, 20 tests	0.2 - 3 ppb

All Glyphosate and Algal Toxins tests are available in a 5-strip format except from item 520046.

RapidScan ST5-W Lateral Flow Assay Reader



Instrument Specifications

Specifications	
Product Name	RapidScan ST5-W Lateral Flow Assay Reader
Article No.	LFRSCAN002
Supplied Materials	RapidScan ST5-W Reader Power Cable and International Adapters QR Code Adapter Check Chart User Guide
Item Dimensions	11 x 14 x 21cm
Item Weight	930 g
Display	5" LCD touch screen
Platform	Lateral Flow
Power	DC 12V 1.5A
Interface	USB 2.0 port
Application Adapters	See table hereunder. Reader price includes one application adapter of choice.

Application Adapters

Article No.	Product Name
LFA0020001	RapidScan ST5-W LFD Reader Adapter for pesticide (glyphosate) and algal toxin strips
LFA0030001	RapidScan ST5-W LFD Reader Adapter for allergen strips
LFA0060001	RapidScan ST5-W LFD Reader Adapter for GMO strips